Salvia apiana In Vitro Shoot System as a source of unique volatile fraction

Agata Krol¹, Adam Kokotkiewicz¹, Bozena Zabiegala², Maria Luczkiewicz¹

¹Department of Pharmacognosy, Faculty of Pharmacy, Medical University of Gdańsk, Gdańsk, POLAND

²Department of Analytical Chemistry, Faculty of Chemistry, Gdańsk University of Technology, Gdańsk, POLAND



Introduction

White sage (Salvia apiana Jepson) is an endemic species typical for the chaparral plant formation of mild, Mediterranean-type climate area of the California Floristic Province in North America. S. apiana has been used as a traditional medicinal and ritual plant by the Native North American Chumash people for a long time. Aerial parts of white sage have been used in folk medicine to treat inflammation and as a calmative agent (Krol et al., 2021). The aim of the study was to established and scale-up S. apiana microshoot cultures in order to obtain sustainable and stable source of S. apiana essential oil (EO).

Plant Materials and Methods



In vitro culture initiation

Seeds sterilization (5.25% NaClO 10 min)



In vitro culture optimization

Growth profile and growth parameters determination



Scaling up microshoots

Cultivation in RITA, Plantform and spray glass bioreactor (SGB)



Essential oil isolation and analysis

The isolation of *S. apiana in vitro* biomass essential oil in Clevenger apparatus and GC-FID/MS analysis. Comparative analysis of the *S. apiana* and *S. officinalis* raw material (RM)

Conclusions

- For the first time, S. apiana in vitro shoot cultures were established as the alternative, continuous source of volatile fraction.
- The largest biomass accumulation and essential oil content were achieved with application of prototype spray glass bioreactor for 4 weeks.
- The essential oil was rich in 1,8-cineole, α -pinene, β -pinene and did not contain thujone.

Results




Gi(growth index) = 641% EO content = 1.27% Schenk-Hildebrandt medium with Thidiazuron 0.22 mg/l and 6benzylaminopurine (2.0



RITA bioreactor

28 days Gi = 641% EO content = 1.10%



Growth profile of S. apiana microshoots, cultivated in RITA® bioreactor, in comparison with the agitated liquid culture



Plantform bioreactor


28 days Gi = 521% EO content = 1.17%



Spray glass bioreactor

28 days Gi = 576%

EO content = 1.27%



References

Krol, A., Kokotkiewicz, A., & Luczkiewicz, M. (2021). White Sage (Salvia apiana)—a Ritual and Medicinal Plant of the Chaparral: Plant Characteristics in Comparison with Other Salvia Species. Planta Medica. doi: 10.1055/a-1453-0964
The study was supported by the project POWR.03.02.00-00-I014/17-00 co-financed by the European Union through the European Social Fund under the Operational Programme Knowledge Education Development 2014—2020